



7-SIGMA Simulation Systems

## CRICOID SKIN STRIPS

**The Cricoid Skin Strip** offers accurate representations of the human anatomy for use in cricothyrotomy training. With an emphasis on realism, our skin offers unique tissue layers consisting of the epidermis, hypodermis, and cricothyroid membrane.



### Features & Solutions

- Works with a variety of simulation trainers
- Caucasian and African American skin colors
- Cricothyroid membrane base layer that provides realistic tactile feedback during cutting
- Product dimensions: 18" x 2.50" x 0.2" per strip
- Weight: Approximately 5 oz. per strip

**7-SIGMA Simulation Systems**  
2843 26th Ave S, Minneapolis, MN 55406  
Phone: 888.722.7224 • Fax: 844.898.4755  
Email: customerservice@7-sigma.com



@7S3USA



@7SIGMAMPLS

<https://www.7-S3.com>