



7-SIGMA Simulation Systems

Cricoid Skin Strips

The Cricoid Skin Strip offers accurate representations of the human anatomy for use in cricothyrotomy training. With an emphasis on realism, our skin offers unique tissue layers consisting of the epidermis, hypodermis, and cricothyroid membrane.

Our skin strips are designed to work with a variety of simulation trainers and offers an 18" x 2.75" cutting area.



Features & Solutions

- High fidelity skin colored epidermis layer
- Caucasian and African American skin colors
- Hypodermis layer complete with fatty tissue and blood vessels
- Cricothyroid membrane base layer that provides realistic tactile feedback during cutting
- Product dimensions are 18" x 2.75" x 0.2" per strip
- Weight is 5oz per strip

WWW.7-SIGMA.COM/7-S3